Impulsion

BODY MASSAGER

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We would like to thank you for purchasing a BODY MASSAGER by IMPULSION. BODY MASSAGER is an appliance that is ideal for making you feel better in your skin. This multifunction massaging device not only helps you prevent unsightly fat deposits but also to reduce the «orange peel» effect. It also contributes to improving blood circulation and softening your skin. In addition, BODY MASSAGER may also be used as a body massaging system.

PLEASE READ ALL THE INSTRUCTIONS PRIOR TO USING YOUR «BODY MASSAGER», IN PARTICULAR THESE FEW BASIC SAFETY INSTRUCTIONS:

- The BODY MASSAGER should only be used for the use for which it is intended and which is described in this manual.
- Your massaging device is designed exclusively for home use.
- In order to ensure that the device continues to operate correctly, we recommend using BODY MASSAGER for 20 minutes and leaving a gap of 10 minutes between each massaging session.
- Please ensure that you do not use the massaging device on the following parts of your body:
 - any part of your body that is swollen, burnt or inflamed or on any part of your body that has a rash or ulcers.
 - A painful abdomen.
 - Legs suffering from varicose veins.
 - Sensitive genital areas.
- People suffering from diabetes and pregnant women should consult their doctor before using the BODY MASSAGER.
- The sensation achieved by the massage should be pleasant and comfortable. If it becomes painful or uncomfortable, please stop all use and consult a doctor.
- This appliance is not intended for use by persons (including children) whose physical, sensorial or mental capabilities are impaired, or by persons lacking experience or knowledge, unless they are supervised or instructed in the use of the appliance first by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
- Make sure you do not drop the device in order to avoid damaging its internal components.
- Never submerge the appliance in water, in your bath or under the shower.
- In order to avoid overheating, refrain from using the equipment under a blanket or pillow.
- If you are not going to use this appliance for a long period, please remove the batteries in order to avoid oxidation.
- If the batteries have become corroded, replace them and clean the battery compartment with a dry cloth.

FUNCTIONS AND USE

1. Rollers
2. Operating light
3. On-off switch
4. Cover for batteries (not supplied)

Use your BODY MASSAGER daily to prevent fat deposits and reduce your skin's «orange peel effect». You will obtain the best results by combining this treatment with a balanced diet and physical exercise.

The BODY MASSAGER has been designed especially for women who take care of their body.

In fact, many women suffer from problems linked to fat deposits and water retention, which leads to the appearance of cellulite which gives their skin this «orange peel» look.

The «orange peel» look is generally used to describe areas of the body where the skin is bumpy (thighs, buttocks, upper arms, etc.). Sometimes cellulite is not visible to the naked eye and is only visible when the skin is pressed between your fingers.

INSTRUCTIONS BEFORE USE

CHANGING THE BATTERIES

- Open the battery cover (4) and insert 2 AA 1.5V batteries (not supplied with the appliance) making sure you keep to the polarity. Then replace the battery cover.
- Turn the switch (3) to the «On» position. The red indicator (2) comes on.

NOTE: No period of use should exceed 20 minutes. Leave the device switched off for 10 minutes after each session.

HOW TO USE YOUR MASSAGING DEVICE

This massaging system is equipped with two soft rollers that help combat cellulite. This system is equipped with an internal motor that generates a vibrating massage function for your muscles and your skin.

We recommend regular use of the BODY MASSAGER in order to achieve the best results:

- 5 minutes daily for each area that requires treatment with a total period of 20 minutes per session.
- You should always massage towards your heart.

NOTE:

Your skin may be slightly red and warm after treatment (this indicates that blood circulation has been stimulated).

The appearance of these symptoms is normal after each massage session and they will disappear quickly. If you wish you can apply a soothing lotion (not supplied) after treatment.

The results may be different from one woman to another due to differences in terms of cellulite, skin type, lifestyle and habits.

Thanks to the vibrating massage combined with the effect of skin contraction, the BODY MASSAGER has a real effect on cellulite.

MAINTENANCE AND CLEANING

1. Always remove the batteries before cleaning the appliance.
2. Never use liquid or an abrasive cleaner to clean this device. You just need to wipe it with a damp cloth.
3. Clean the battery compartment using a dry cloth.

ADVICE REGARDING THE DISPOSAL OF WASTE:

- All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.
- When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law. Prior to disposal, please remove the battery and dispose of it at a collection point so that it can be recycled.

Please note: under no circumstances should spent batteries be disposed of with household waste!!!

LIMITED WARRANTY

IMPULSION® guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.

The IMPULSION® guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this IMPULSION® product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside IMPULSION®'s control.

IMPULSION® may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, IMPULSION® will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through IMPULSION®'s Service Centre®. Any attempt to maintain this product by a person other than IMPULSION®'s Service Centre will render this guarantee void.

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